

# Pelican Pool & Spa Hoist Solid Seat Instructions

Ensure locking pin is engaged



## BEFORE MOVING HOIST TO POOL DECK (MOBILE VERSION)

1. Install the battery. On receipt of your hoist the battery will be charged and ready to go.
2. **NEVER** remove the hand control or the control box from the hoist.
3. For models with a battery cover, leave on the control box when not in use.

## INSTALLATION AT POOL DECK (MOBILE VERSION)

1. Wheel the hoist to the pool deck and place the base plate of the hoist over the anchor points which have been inserted into the pool deck.
2. Insert the Eye Bolts into the holes and hand tighten only.

## OPERATING THE POOL HOIST

1. **Always** test the operation of the hoist by using the hand control, test the up and down and side to side operation before commencing a lift with a client.
2. Lift the arm rest of the seat in a horizontal motion. Once the client has transferred onto the seat commence the lift.
3. Using the hand control, raise the lifting arm high enough for the client's feet to be off the ground.

4. Using the hand control, manoeuvre the chair over the water, to the landing zone and lower.
5. **NEVER** push the hoist via the actuator to swing over the pool/spa, this will damage the hoist.

## BATTERY CARE

1. A battery charger is supplied and should be placed on a wall near a 240v power point in a weather / vandal proof area.
2. After use, the battery should be placed on the charger to ensure it is fully charged when required next.

## SAFETY AND CARE

1. Push the red button to activate the emergency stop. To release the emergency stop take off the battery and replace it back into the unit again.
2. In case of an emergency activate the manual lowering function of the actuator (see user guide). Rotate the grey or red knob on the top of the actuator clockwise to lower the arm of the hoist.

## TROUBLE SHOOTING

### Hoist won't work?

- Check the emergency stop button is not activated.
- Check the battery is charged.
- Check the battery is installed correctly.
- Check the hand control hasn't been damaged.
- Check for any corrosion on the terminals of the battery and the control box.
- Call Para Mobility Pelican® Installation & Service or your supplier if you have any questions or concerns.



**Para Mobility**

Making everyday living easier

**Phone 1300 444 600**

Pelican® Installation & Service  
Delivery, installation + servicing Australia wide  
enquiries@paramobility.com.au  
www.paramobility.com.au