



# PLATYPUS MOBILE CHANGE TABLE



This change table as a lightweight alternative to standard change tables. The unique design allows a person to be changed in any area (wet or dry) whilst being completely mobile.

The Platypus Mobile Change Tables has two different designs – a flat bed, or an adjustable back.

## FEATURES

- Sturdy frame
- 700mm (wide) x 1700mm (long) x 800mm (high)
- Stainless steel castors (2 lockable)
- Safety belt/s
- Easy to dry and clean
- SWL 80kg
- Accidental burns are prevented as plastic does not absorb the heat
- Open weave fabric which allows the easy passage of air and water
- Rounded surfaces which improve comfort and reduces the risk of skin tears

**Para Mobility**  
*Making everyday living easier since 1983*

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide  
[sales@paramobility.com.au](mailto:sales@paramobility.com.au)

[www.paramobility.com.au](http://www.paramobility.com.au)

