



KINGFISHER™ POOL HOIST



Independence and access to the water! Compact and adaptable.

The Kingfisher™ Pool Hoist, designed and manufactured in Australia by Para Mobility, will allow access to a pool where traditional methods are not able to be used. Designed in close consultation with Occupational Therapists and other industry professionals, the Kingfisher™ Pool Hoist is a simple and safe operation that is both comfortable and dignified.

The Kingfisher™ is ideal for home pools, retirement villages and hotels, plus more. It is user-operated and secured to the pool deck easily. Perfect for people with upper body control and those who can self transfer.

The Kingfisher™ Pool Hoist has a Safe Working Load (SWL) of 150kg, rechargeable Linak 24v battery, floating hand control, load bearing lift-up arm rests for ease of transfer, removable and adjustable foot rest, flip up seat for storage, wheel away transport trolley option for storage, lap belt for safety, left or right side mounting, extra heavy duty Linak Actuator with emergency lifting.



KINGFISHER™ POOL HOIST **PRODUCT SNAPSHOT**

SPECIFICATIONS

- **Lifting Capacity:** SWL 150kg
- **Weight of Hoist:** 48kg (without the battery)
- **Max—Deck to water level:** Up to 250mm
- **Max—Set Back:** 500mm
- **Stow away size:** H 1230 x W 550 x D 400
- **Seat size:** 400mm x 400mm
- **Back rest size:** 400mm 300mm
- **Foot rest size:** 300mm x 150mm (can be modified)
- **Powder coated:** White

STANDARD INCLUSIONS

- Lift up arm rests
- 24v Battery
- Battery charger
- Battery cover
- Hand control
- Seat belt
- All weather cover
- Mounting kit

OPTIONAL ACCESSORIES

- Foot rest
- Chest strap
- Transport trolley

MATERIALS

- **Base Plate:** 316g Stainless Steel
- **Lifting Frame:** 316g & 304g Stainless Steel
- **Actuator:** Linak
- **Battery:** Linak
- **Control Box:** Linak
- **Floating Hand Control:** Linak
- **Seat:** High Density Polyethylene
- **Seat Belt:** Interwoven nylon/polyester
- **Foot rest:** High Density Polyethylene



View in action:

www.paramobility.com.au/mobility-products/outdoor-products/kingfisher-poolspa-hoist

Para Mobility
Making everyday living easier since 1983

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide
sales@paramobility.com.au

www.paramobility.com.au



Registered NDIS Provider

