



FREE STANDING TRACK (FST)

When safety comes first

The Free Standing Track provides a more effective transfer system than the traditional mobile lifter. It is lightweight and economical, providing the utmost in safety for both client and caregiver.

Ideal for travelling and great for rental properties where a permanent ceiling mount hoist is not an option. The Free Standing Track offers the freedom to move from room to room.



FREE STANDING TRACK (FST) PRODUCT SNAPSHOT

HIGHLIGHTS

- Total system weighs approximately 20kg making it the lightest system on the market
- No single component weighs more than 7.5kg
- Safe working load of 136kg
- 2.4m length allows ample room for transfers from bed to chair
- Secure fit supports ensure stability
- Assembly in minutes by one person with no tools required
- Components fold easily to fit in the ergonomic carry bag

SPECIFICATIONS

- Construction: Aluminium
- Height Adjustable: up to 1.95 m
- Length: 2.4 m
- Unit Weight: 20kg total
- Maximum Load: 136 kgs
- Bag dimensions: L1480 x W 410 x H320MM
- Tested to ISO standards

NOTE: Hoist carriage and lifting arm not included.


Para Mobility
Making everyday living easier since 1983

Phone (02) 9651 4446

Delivery, installation + servicing Australia wide
sales@paramobility.com.au

www.paramobility.com.au

