



MOBILE HEIGHT ADJUSTABLE CHANGE TABLE



The freedom to safely wash and change

Designed and manufactured in Australia by Para Mobility, the Mobile Height Adjustable Change Table allows easy transfer from a wheelchair when the table top is at its lowest point.

The table can then be raised to a working height for the carer. The table top is removable and can be tailor made to suit the dimensions of the room and the user.

NB: Care should be taken in a wet environment.



MOBILE HEIGHT ADJUSTABLE CHANGE TABLE

PRODUCT SNAPSHOT

FEATURES

- Lowers to allow a person to slide transfer or use a floor lifter from a wheelchair to the table top
- Hand held control for easy operation
- Table top can be easily removed for easy storage, transport or fitting into tight areas
- Top is padded and covered with high grade PVC material for urine and stain control
- Unit is on castors with 2 locking wheels
- Framework is galvanised and powder coated for corrosion resistance
- Used for washing and changing in either a wet or dry area, although not suitable for use under a shower

SPECIFICATIONS

- Restraint strap included
- Stainless steel shelf
- Standard table top size: 1850mm long x 700mm wide. Sizes can be made to suit the user.
- Rises to approximately 880mm
- Lowest point approximately 480mm
- Has clearance of 150mm at bottom allowing floor hoist to fit underneath
- SWL (safe working load) 150kg

OPTIONAL EXTRAS

- Drop down sides 150mm standard height x 1200 long. Sides can be made to order.
- Battery power pack, allowing for use without being plugged into a power point
- Additional straps
- Foot Pedal operation



Phone (02) 9651 4446

Delivery, installation + servicing Australia wide
sales@paramobility.com.au

www.paramobility.com.au



Registered NDIS Provider

