



# IBIS® Portable Car Access Lifter



*Freedom to explore*

The IBIS® Portable Car Access Lifter system has been designed and manufactured in Australia by Para Mobility to allow people with limited mobility the option of getting in a car, SUV, 4WD or van. Passengers can be seated in the front passenger seat of most cars, 4WD vehicles and the side entry of vans and caravans. It is operated by the person or carer, with a press button action allowing the person to be transferred from a wheelchair into the vehicle without any lifting, whilst being completely portable.

To further make the IBIS® Access Lifter even more portable an IBIS® Car Bracket is available, made to attach to the hinges of your car door, which eliminates the need for any permanent car modifications. Clip the IBIS® Car Bracket onto the door hinge of the car, attach the IBIS® Portable Car Access Lifter to the bracket as normal, plug in the power source using the portable battery pack and you are ready to go.

## FEATURES

- TGA approved
- NDIS confirmed
- Operation can be by either the user or the carer
- Press button hand control which is powered from the portable battery pack or vehicle battery
- Easily removed and stored in the boot when not in use and whilst travelling
- Lightweight for easy handling – 11kg
- Includes IBIS® sling – small, standard or high back (tailor made slings can also be arranged)
- 3 x size options suitable for small, medium and large vehicles

