


Kingfisher® Pool Access Seat



The Kingfisher® Pool Access Seat, designed and manufactured in Australia by Para Mobility, will allow access to a pool where traditional methods are not able to be used.

Designed in close consultation with Occupational Therapists and other industry professionals, the Kingfisher® Pool Access Seat is a simple and safe operation that is both comfortable and dignified.

The Kingfisher® is ideal for home pools, retirement villages and hotels, plus more. It is user-operated and secured to the pool deck easily. Perfect for people with upper body control and those who can self transfer.



FEATURES

- TGA approved
- NDIS confirmed
- Lift up arm rests
- 24v Battery
- Battery charger
- Battery cover
- Hand control
- Seat belt
- All weather cover
- Transport trolley
- Mounting kit

SPECIFICATIONS

- Lifting Capacity: SWL 150kg
- Weight of Hoist: 48kg (without the battery)
- Max—Set Back: 500mm
- Stow away size: H 1230 x W 550 x D 400
- Seat size: 400mm x 400mm
- Back rest size: 400mm 300mm
- Foot rest size: 300mm x 150mm (can be modified)
- Powder coated: White
- Base Plate: Stainless Steel
- Lifting Frame: Stainless Steel
- Linak actuator, battery, control box and floating hand control
- High Density Polyethylene seat and foot rest
- Manufactured to Australian Standard — AS3581 and eye bolts

OPTIONAL EXTRAS

- Foot rest
- Chest strap

Freedom to explore

