

- Weight: 6.7kg w/o spreader bar
- Measurements: 320 x 130mm
- Battery operation: 36V lead acid batteries
- Charging: 100–230V. 0.5A
- Chassis: Aluminium
- Housing: High impact, moulded ABS
- Noise: LWA during lifting/lowering without load is 49dB (A)
  LWA during lifting/lowering at maximum load is 53dB (A)
- Lifting speed down: 44mm/ sec with 89kg load
- · Lifting speed up: 30mm/ sec with 89kg load
- No. of lifts per charge: 50 lifts of 0.5m with 100kg
- ISO: ISO 10535



## Phone 1300 444 600

Delivery, installation + servicing Australia wide enquiries@paramobility.com.au www.paramobility.com.au

**LEARN MORE** 

