


PELICAN POOL HOIST MEASUREMENT GUIDE

CUSTOMER NAME	CUSTOMER PHONE
POOL HOIST COLOUR	
NOTES	
	TO YOUR POOL AND MEASURE WHERE APPLICABLE.
DO YOU HAVE THE MINIMUM CONCRETE WIDTH AN	ID THICKNESS TO MOUNT THE HOIST ON POOL DECK? Y
Pelican Rim Setback (Aerial View) Pool/Spa Edging Pelican Setback Pelican Setback 100 mm	WATER DEPTH - IMPORTANT NOTICE In accordance with the Disabilities Standards Code, the Pelican base rim must be positioned on the pool or spa surrounds where the water dept is not less than 500mm and not more than 1300mm.
Pelican setback is the standard setback distance (100mm).	
Pool Depth at Location (Example 1.2m deep should be expresse	ed as 1200)
Round measurer	ments to the nearest 10mm
Profile A 1 Coping	
To Water Level Only 2 Water Level Pool Deck	A2 INTERNAL RISER A3 EXTERNAL RISER
	NOTE: Measurement 2 may be different to measurement 3.
Profile B To Water Level Only Water Level	B2 INTERNAL RISER
Profile C 1 Grate Pool Deck	C1 GRATE WIDTH C2 INTERNAL RISER
Water Level	





PELICAN POOL HOIST COLOUR SELECTION GUIDE

PLEASE SELECT THE COLOUR YOU WOULD LIKE YOU PELICAN POOL HOIST TO BE.



