



# Aquatic Pool Wheelchairs



The Aquatic Pool and Wet Area Wheelchair, designed and manufactured in Australia by Para Mobility, is a lightweight, rustproof, self-propelled wheelchair designed to be fully immersed in salt, chlorinated or fresh water. This wheelchair is made from quality furniture grade PVC, which will not absorb the heat if left out in the sun, while the strong open weave fabric allows passage of air and water. The chair has solid tyres and is ideal for leisure centres, hydrotherapy pools, showering and outdoor environments. The Platypus and Pelican Pool and Wet Area Wheelchairs have been designed for comfort, safety and practicality, and are easily and safely self-propelled. The wheels are equipped with hand brakes on each side of the wheelchair for maximum safety and control. The brakes are strategically positioned to allow carers easy access, at hip height, when activating them. An option is available for the aquatic wheelchair to be user controlled giving them direct access and independence.

The Platypus model is used for on ground and ramps while the Pelican model includes a lifting frame and lifting straps to be used in conjunction with Pelican Pool Hoist.

## FEATURES

- NDIS confirmed
- Robust fixed plastic frame
- 600mm rear wheels
- Hand brakes (with upgrade option for dual control)
- Swing away arms for side transfer
- High quality stainless steel castors
- Folding footrest for front entry
- Strong PVC construction (does not get hot)
- The strong open weave fabric allows a flow of air or water
- Easy to dry and clean
- Lightweight frame avoiding carer injuries



*Freedom to explore*



**Para Mobility**  
Making everyday living easier

**Phone 1300 444 600**

Delivery, installation + servicing Australia wide  
[enquiries@paramobility.com.au](mailto:enquiries@paramobility.com.au)  
[www.paramobility.com.au](http://www.paramobility.com.au)

LEARN MORE

