

Pelican Plus[®] Pool & Spa Hoist



Building on the trusted design of the standard Pelican[®] Pool and Spa Hoist, the Pelican Plus[®] Pool and Spa Hoist incorporates an advanced full 360-degree powered rotation motor, allowing users to control both the up-and-down lifts and the left-and-right pivot motions with ease.

Designed and manufactured in Australia by Para Mobility, the Pelican Plus is the next evolution in pool access solutions, offering unparalleled independence and autonomy to users with disabilities and mobility limitations.

FEATURES

- TGA approved and NDIS confirmed
- Easily wheeled from one area to another for multiple use and or storage (mobile unit)
- No manual lifting or rotating required
- Powered rotation motor enables seamless 360-degree movement controlled via the four direction hand control
- Standard lifting capacity 150kgs and 200kg
- Cover included
- Easily installed
- Complies with Disability Discrimination Act (Access to Premises – Buildings) under the Building Code of Australia

SEATING OPTIONS

- Floating sling
- Solid plastic seat
- Pelican Aquatic Wheelchair – designed and manufactured by Para Mobility, a great option as it can be used around the pool area
- Platform lift (Meets Access to Premises Standards 2010 for pools greater than 70m perimeter)

OPTIONAL EXTRAS

- Earthing ability
- Bluetooth handset

SPECIFICATIONS

- Safe Working Load: 150kg and 200kg
- Batteries: Linak – 24v re-chargeable battery
- Rotation: 360 degree power rotation
- Total Weight: From 60kg
- Floor Mounting Plate: 450 diameter
- Construction: Stainless Steel, powder coated to colour of choice
- Hand control: Linak – press-button four directional control hand switch
- Installation: Anchored to pool deck with Stainless Steel anchors and eye bolts



Freedom to explore

