



Para Mobility

Making everyday living easier

 **horizon.**



Montar Wall Mounted Change Table

Product Code: 7-MONTAR

The Montar Wall Mounted Change Table offers space saving efficiencies and will integrate seamlessly into your facility. It includes state of the art electrics to achieve a 385mm low height and an 985mm high height for ergonomic transfer of the patient.


- Charcoal Vinyl Padded Top
- Corded Handset
- Drop Down Safety Rail Folds Up For Saving Space
- Electric Hi-Low
- Wall Mounted
- APO Grey Powdercoating
- Battery Back Up

Optional Features - Added at time of manufacture

- Custom Width & Length
- Plas Coat Coating Upgrade
- Bluetooth Handset
- Custom Vinyl Colour
- Custom Powdercoat Colour

Accessories: (Added after manufacture)

- Lap belt



Specifications:

| Size | Paediatric | Adult |
|-------------------------------|---------------|---------------|
| Sleeping Surface | 700 x 1400 mm | 700 x 1880 mm |
| Overall Length | 1400 mm | 1880 mm |
| Overall Width (Open / Closed) | 350 / 1030 mm | 350 / 1030 mm |
| Height Range (Approximate) | 385 - 985 mm | 385 - 985 mm |
| SWL | 200 kg | 200 kg |